

## Supplementary File 3: Nuclei Area Quantification

for “Four-Fold Expansion Microscopy of Cultured Cells” by Corban Swain, Sapna Sinha, and Ed Boyden

*Jupyter Notebook (exported to .pdf) containing ImageJ commands and Python script for quantification of nuclei area from images.*

### 1. Setting Up Repository and Opening the Notebook

1. You will need to install and set up the following software before continuing:

- ☐ **git** <https://git-scm.com/>
- ☐ **mamba** <https://mamba.readthedocs.io>
- ☐ **Fiji (ImageJ)** <https://imagej.net/software/fiji/downloads> (version 1.54p was used)

2. Clone the repository at <https://github.com/CorbanSwain/ExM-Demo-Figures> into a directory (e.g. ~/ExM-Demo-Figures) on your computer. Then initialize and update the submodules.

```
git clone https://github.com/CorbanSwain/ExM-Demo-Figures ~/ExM-Demo-Figures
git submodule update --init
```

3. Set up the Python environment according to the environment configuration provided in the repository.

```
cd ~/ExM-Demo-Figures/envs/
mamba create -f exm-demo-fig-py3.12-env.yml
```

4. Then run JupyterLab and open this notebook via the JupyterLab interface.

```
mv ~/supplementary_file_01-nuclei-area-quant.ipynb ~/ExM-Demo-Figures/source/
cd ~/ExM-Demo-Figures/source
mamba run -n exm-demo-fig-py3.12-env jupyter lab supplementary_file_03-nuclei-area-quant.ipynb
```

### 2. Segmentation and Measurement of Images in Fiji

Follow these instructions to process the images for nuclei area quantification.

#### 1. Open Image

- Run Fiji and open the image located at ~/ExM-Demo-Figures/data/nuc-size-analysis/expanded-hex-cells-a.tif  
*This is fluorescent microscope image of expanded HEK cells taken with a 40x objective on a Nikon CSU-W1. The DAPI channel has already been separated as its own image file.*
- Image > Adjust > Brightness/Contrast...
  - Click Reset
  - Click Auto
- LUT Toolbar Button > Grays

#### 2. Duplicate Image


- Image > Duplicate...
- Title: = expanded-hek-cells-a-MASK.tif
- Perform the following steps on the expanded-hek-cells-a-MASK.tif file.

### 3. Histogram Normalization

- Process > Enhance Local Contrast (CLAHE)
- blocksize: = 255 for expanded samples (63 for unexpanded)
- histogram bins: = 256
- maximum slope: = 2.00
- mask: = \*None\*
- fast = NO

### 4. Thresholding

- Image > Adjust > Threshold...
- lower threshold value: Adjust until approximately 25% of pixels are selected. *This target percentage may need to be adjusted based on cell confluency.*



- upper threshold value: maximum
- method: Default
- color: Red
- Dark background = YES
- Don't reset range = YES
- 16-bit histogram = NO
- Stack histogram = NO
- Raw values = NO
- Finally, click Apply.

### 5. Binary Closing

- Process > Binary > Options...

- Iterations: = 4 for expanded samples (1 for unexpanded)
  - Count: = 1
  - Black background = YES
  - Pad edges when eroding = YES
  - EDM output: = 8-bit
  - Do: = Close
6. **Binary Fill Holes**
- Process > Binary > Fill Holes
7. **Watershed Segmentation**
- Process > Binary Watershed
8. **Save Mask**
- File > Save...
  - Save in: ~/ExM-Demo-Figures/data/nuc-size-analysis/
  - File name: expanded-hek-cells-a-MASK.tif
9. **Analyze**
- Analyze > Analyze Set Measurements...
    - Area = YES
    - *all others* = NO
    - Decimal Places: = 3
  - Analyze > Analyze Particles...
    - Size (micron<sup>2</sup>): = 500-4000 for expanded samples (31-250 for unexpanded)
    - Pixel units = NO
    - Circularity: = 0.5-1.00
    - Show: = Overlay
    - Display results = YES
    - Clear results = YES
    - Summarize = NO
    - Add to Manager = NO
    - Record starts = NO
    - Exclude on edges = YES
    - Include holes = YES
    - Overlay = YES
    - Composite ROIs = NO
10. **Save Results**
- Select the window Results
  - File > Save As (Within the Results window, not the Fiji toolbar)
  - Save in: ~/ExM-Demo-Figures/data/nuc-size-analysis/
  - File name: expanded-hek-cells-a-RESULTS.csv
11. **Transfer Regions**
- Select the window expanded-hek-cells-a-MASK.tif
  - Image > Overlay > To ROI Manager *the ROI Manager window should now be open.*

- Select the window `expanded-hek-cells-a.tif`
- Select the window **ROI Manager**
- Click the first region (i.e. 1-1) in the ROI manager region list.
- Click **Update**. *All of the regions should now be overlaid on the original image.*

## 12. Inspect Regions

- Use the ROI Manager to view each region in the original image by selecting each region label from the ROI Manager's list or by clicking the region in the image. If a region does not appear to be a complete, non-dividing nucleus, select it and click **Delete**, then continue to the next region.
- Once all the regions have been reviewed in this way, transfer the indices of the remaining good regions into the Python list below (section 3.3.3), setting the variable named `good_roi_exm_a_list`. The expected result for `expanded-hek-cells-a` is

```
good_roi_exm_a_list = range(2, 28) # regions 2 - 27 are good; region 1 is a partial nucleus
```

## 13. Annotate Original Image

- In the ROI Manager, click **Deselect**
- Click **Properties...**
  - **Range:** = *should automatically be filled in with the range of all the regions (e.g 1-26)*
  - **Position:** = none
  - **Group:** = none
  - **Stroke color:** = yellow
  - **Width:** = 4
  - **Fill color:** = none
  - **List coordinates** = NO
- In the main Fiji toolbar, select **Image > Overlay > Labels...**
  - **Color:** = yellow
  - **Font size** = 24 *adjust as needed*
  - **Show labels** = YES
  - **Use names as labels** = NO
  - **Draw backgrounds** = NO
  - **Bold** = NO
- In the ROI Manager, click **Deselect**
- In the ROI Manager, click **Flatten [F]** *this should open a new image window*
- In the main Fiji toolbar, select **File > Save As > PNG...** to save the new image
  - **Save in:** `~/ExM-Demo-Figures/data/nuc-size-analysis/`
  - **File name:** `expanded-hek-cells-a-ANNOTATED.png`

Repeat this process for each of the other nuclei images to be analyzed: `expanded-hek-cells-b.tif` and `normal-hek-cells.tif`.

## 3. Analysis, Plotting of Nuclei Sizes with Python

Run the following cells to process the extracted results data.

### 3.1 Package Imports

```
# 3.1.1 base and metaprogramming imports
import os
import sys
import pandas as pd
import numpy as np
import seaborn as sns
import matplotlib as mpl
import matplotlib.pyplot as plt

# 3.1.2 ipynb-specific imports
try:
    import ipynbname
except Exception as e:
    warnings.warn(f'`ipynbname` module import failed with error: {e!r}')
    ipynbname = None

if ipynbname is not None:
    _file = ipynbname.path()
else:
    # manually set the full path to the current notebook below if `ipynbname`
    # import fails
    _file = ''

# 3.1.3 key path variables and management
source_dir = os.path.split(_file)[0]
repo_dir = os.path.abspath(os.path.join(source_dir, '..'))
data_dir = os.path.join(repo_dir, 'data')
figure_dir = os.path.join(repo_dir, 'figures')
if source_dir not in sys.path:
    sys.path.append(source_dir)
```

### 3.2 Plotting Defaults

```
sns.set_style(
    'ticks',
    {'font.sans-serif': [
        'Helvetica',
        'Arial'],
    'mathtext.fontset': ['cm'],,
```

```

    'axes.linewidth': 2})

axes_lw = 1.25
mpl.rcParams['axes.linewidth'] = axes_lw
mpl.rcParams['xtick.major.width'] = axes_lw
mpl.rcParams['xtick.minor.width'] = axes_lw
mpl.rcParams['ytick.major.width'] = axes_lw
mpl.rcParams['ytick.minor.width'] = axes_lw

mpl.rcParams['axes.facecolor'] = 'none'

mpl.rcParams['figure.frameon'] = False
mpl.rcParams['figure.facecolor'] = 'none'
mpl.rcParams['figure.edgecolor'] = 'none'

```

### 3.3 Filepaths and Other Parameter Inputs

```

# 3.3.1 filepath names
# Edit these fields to match the files you wish to analyze.
nuc_size_dname = 'nuc-size-analysis'
normal_fname = 'normal-hek-cells-RESULTS.csv'
exm_a_fname = 'expanded-hek-cells-a-RESULTS.csv'
exm_b_fname = 'expanded-hek-cells-b-RESULTS.csv'

figure_name = 'nuc_size_comparison' # output figure name

# 3.3.2 table keys
blank_index_key = ' ' # Fiji populates the results `.csv` with a blank
                     # header for the first/index column
index_key = 'original_index'
category_key = 'category'
expanded_cat_name = 'expanded'
unexpanded_cat_name = 'normal'
file_key = 'original_file'
area_key = 'Area'
sqrt_area_key = '[Area]^(1/2)'

# 3.3.3 good region list (from manual inspection)
good_roi_normal_list = [
    9,  13, 15, 17, 19, 20, 21, 26, 28, 29,
    30, 31, 35, 36, 37, 38, 39, 40, 41, 43,

```

```

44, 45, 46, 47, 48, 49, 51, 54, 55, 57,
58, 59, 60, 61, 63, 65, 68, 69, 70, 75,
76, 78, 79, 81, 82, 83, 84, 85, 86, 87,
90, 92, 94, 96, 99, 103, 104, 105, 106, 109,
111, 112, 113, 114, 115, 116, 117, 118, 119, 120,
124, 128, 134, 137, 138, 142, 143, 144, 145, 148,
150, 152, 153, 154, 155, 156, 163, 164, 166, 170,
173, 174, 176, 177, 178, 180, 181, 184, 186, 188,
189, 190, 191, 192, 193, 194, 198, 199, 200, 201,
202, 204, 207, 209, 212, 213, 214, 215, 219, 220,
221, 222, 224, 225, 228, 231, 234, 235, 236, 237,
241, 242, 243, 246, 247, 249, 252, 254, 260, 264,
265, 267, 268, 270, 280, 281, 287, 288]
good_roi_exm_a_list = range(2, 28) # 2 to 27, inclusive
good_roi_exm_b_list = range(1, 29) # all regions, 1 to 28, inclusive

```

### 3.4 Loading and Processing Data Tables

```

# 3.4.1 generate full paths to files
nuc_size_dir = os.path.join(data_dir, nuc_size_dname)
normal_path = os.path.join(nuc_size_dir, normal_fname)
exm_a_path = os.path.join(nuc_size_dir, exm_a_fname)
exm_b_path = os.path.join(nuc_size_dir, exm_b_fname)

# 3.4.2 load `.csv` files into pandas `DataFrame`s
raw_normal_df = pd.read_csv(normal_path)
raw_exm_a_df = pd.read_csv(exm_a_path)
raw_exm_b_df = pd.read_csv(exm_b_path)

# 3.4.3 process the raw tables
def process_raw_df(raw_df, fname, cat_name, good_roi_list):
    df = raw_df.copy()
    df[file_key] = fname
    df[category_key] = cat_name
    df = df.rename(columns={blank_index_key: index_key})
    df = df[df[index_key].isin(list(good_roi_list))]
    return df.reset_index(drop=True)

normal_df = process_raw_df(
    raw_normal_df, normal_fname, unexpanded_cat_name, good_roi_normal_list)
exm_a_df = process_raw_df(

```

```

    raw_exm_a_df, exm_a_fname, expanded_cat_name, good_roi_exm_a_list)
exm_b_df = process_raw_df(
    raw_exm_b_df, exm_b_fname, expanded_cat_name, good_roi_exm_b_list)

# 3.4.4 make a concatenated table of all data points
nuc_size_df = (
    pd.concat([exm_a_df, exm_b_df, normal_df])
    .reset_index(drop=True))
nuc_size_df[category_key] = nuc_size_df[category_key].astype('category')
nuc_size_df[sqrt_area_key] = np.sqrt(nuc_size_df[area_key])

# 3.4.5 log table information and preview
def print_df_info(df, name):
    print(f'\n{name} Data Frame: ' + '-'*70)
    df.info(memory_usage=False)
    print()
    print(df.head(2))

def print_df_infos(raw_df, df, name):
    print_df_info(raw_df, 'Raw ' + name)
    print_df_info(df, name)

print_df_infos(raw_normal_df, normal_df, 'NORMAL')
print_df_infos(raw_exm_a_df, exm_a_df, 'EXPANDED-A')
print_df_infos(raw_exm_b_df, exm_b_df, 'EXPANDED-B')
print_df_info(nuc_size_df, 'ALL')

```

```

Raw NORMAL Data Frame: -----
<class 'pandas.core.frame.DataFrame'>
RangeIndex: 297 entries, 0 to 296
Data columns (total 2 columns):
#   Column  Non-Null Count  Dtype
---  -
0      297 non-null    int64
1   Area    297 non-null    float64
dtypes: float64(1), int64(1)
      Area
0  1  51.732
1  2  68.362

```

```

NORMAL Data Frame: -----

```



```
<class 'pandas.core.frame.DataFrame'>
RangeIndex: 148 entries, 0 to 147
Data columns (total 4 columns):
#   Column          Non-Null Count  Dtype
---  -
0   original_index  148 non-null   int64
1   Area            148 non-null   float64
2   original_file    148 non-null   object
3   category         148 non-null   object
dtypes: float64(1), int64(1), object(2)
      original_index      Area      original_file  category
0              9  68.995  normal-hek-cells-RESULTS.csv  normal
1             13  92.185  normal-hek-cells-RESULTS.csv  normal
```

Raw EXPANDED-A Data Frame: -----

```
<class 'pandas.core.frame.DataFrame'>
RangeIndex: 27 entries, 0 to 26
Data columns (total 2 columns):
#   Column  Non-Null Count  Dtype
---  -
0      27 non-null   int64
1   Area  27 non-null   float64
dtypes: float64(1), int64(1)
      Area
0  1  1377.171
1  2  1163.395
```

EXPANDED-A Data Frame: -----

```
<class 'pandas.core.frame.DataFrame'>
RangeIndex: 26 entries, 0 to 25
Data columns (total 4 columns):
#   Column          Non-Null Count  Dtype
---  -
0   original_index  26 non-null   int64
1   Area            26 non-null   float64
2   original_file    26 non-null   object
3   category         26 non-null   object
dtypes: float64(1), int64(1), object(2)
      original_index      Area      original_file  category
0              2  1163.395  expanded-hek-cells-a-RESULTS.csv  expanded
1              3  1108.642  expanded-hek-cells-a-RESULTS.csv  expanded
```

```
Raw EXPANDED-B Data Frame: -----
<class 'pandas.core.frame.DataFrame'>
RangeIndex: 28 entries, 0 to 27
Data columns (total 2 columns):
#   Column  Non-Null Count  Dtype
---  -
0      28 non-null    int64
1   Area  28 non-null    float64
dtypes: float64(1), int64(1)
      Area
0 1   965.444
1 2  1149.873
```

```
EXPANDED-B Data Frame: -----
<class 'pandas.core.frame.DataFrame'>
RangeIndex: 28 entries, 0 to 27
Data columns (total 4 columns):
#   Column          Non-Null Count  Dtype
---  -
0   original_index  28 non-null    int64
1   Area            28 non-null    float64
2   original_file    28 non-null    object
3   category         28 non-null    object
dtypes: float64(1), int64(1), object(2)
      original_index      Area      original_file  category
0      1      965.444  expanded-hek-cells-b-RESULTS.csv  expanded
1      2  1149.873  expanded-hek-cells-b-RESULTS.csv  expanded
```

```
ALL Data Frame: -----
<class 'pandas.core.frame.DataFrame'>
RangeIndex: 202 entries, 0 to 201
Data columns (total 5 columns):
#   Column          Non-Null Count  Dtype
---  -
0   original_index  202 non-null    int64
1   Area            202 non-null    float64
2   original_file    202 non-null    object
3   category         202 non-null    category
4   [Area]^(1/2)     202 non-null    float64
dtypes: category(1), float64(2), int64(1), object(1)
      original_index      Area      original_file  category  \
0      2  1163.395  expanded-hek-cells-a-RESULTS.csv  expanded
```

```

1          3 1108.642 expanded-hek-cells-a-RESULTS.csv expanded

[Area]^(1/2)
0    34.108577
1    33.296276

```

### 3.5 Show ROI Delete List

This is a tool to print the regions that were removed from the original ROI Manager at step 2.12 above.

```

def print_list_in_columns(data, num_columns=10, column_padding=2):
    # function partially generated by Gemeni AI, prompted by Corban Swain
    max_len = max(len(str(item)) for item in data) if data else 0
    col_width = max_len + column_padding
    for i, item in enumerate(data):
        print(str(item).ljust(col_width), end="")
        if (i + 1) % num_columns == 0:
            print()
    if len(data) % num_columns != 0:
        print() # Print a final newline if the list doesn't perfectly fill the last row

def print_delete_list(good_list, raw_df, name):
    good_len = len(good_list)
    total_len = len(raw_df)
    bad_len = total_len - good_len
    print(f'\nLabels of deleted ROIs for {name} dataset:')
    print(f'({good_len:d} regions kept; '
          f' {bad_len:d} regions deleted)')
    if bad_len == 0:
        print('none')
    return
    print_list_in_columns(
        [(i+1) for i in range(total_len) if (i+1) not in good_list])

print_delete_list(good_roi_normal_list, raw_normal_df, 'NORMAL')
print_delete_list(good_roi_exm_a_list, raw_exm_a_df, 'EXPANDED-A')
print_delete_list(good_roi_exm_b_list, raw_exm_b_df, 'EXPANDED-B')

```

Labels of deleted ROIs for NORMAL dataset:

(148 regions kept; 149 regions deleted)

```

1  2  3  4  5  6  7  8  10  11
12 14 16 18 22 23 24 25 27 32

```

|     |     |     |     |     |     |     |     |     |     |
|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| 33  | 34  | 42  | 50  | 52  | 53  | 56  | 62  | 64  | 66  |
| 67  | 71  | 72  | 73  | 74  | 77  | 80  | 88  | 89  | 91  |
| 93  | 95  | 97  | 98  | 100 | 101 | 102 | 107 | 108 | 110 |
| 121 | 122 | 123 | 125 | 126 | 127 | 129 | 130 | 131 | 132 |
| 133 | 135 | 136 | 139 | 140 | 141 | 146 | 147 | 149 | 151 |
| 157 | 158 | 159 | 160 | 161 | 162 | 165 | 167 | 168 | 169 |
| 171 | 172 | 175 | 179 | 182 | 183 | 185 | 187 | 195 | 196 |
| 197 | 203 | 205 | 206 | 208 | 210 | 211 | 216 | 217 | 218 |
| 223 | 226 | 227 | 229 | 230 | 232 | 233 | 238 | 239 | 240 |
| 244 | 245 | 248 | 250 | 251 | 253 | 255 | 256 | 257 | 258 |
| 259 | 261 | 262 | 263 | 266 | 269 | 271 | 272 | 273 | 274 |
| 275 | 276 | 277 | 278 | 279 | 282 | 283 | 284 | 285 | 286 |
| 289 | 290 | 291 | 292 | 293 | 294 | 295 | 296 | 297 |     |

Labels of deleted ROIs for EXPANDED-A dataset:  
 (26 regions kept; 1 regions deleted)  
 1

Labels of deleted ROIs for EXPANDED-B dataset:  
 (28 regions kept; 0 regions deleted)  
 none

### 3.6 Compute Median Nuclei Sizes and $n$

```

pivot = nuc_size_df.pivot_table(
    columns=category_key,
    values=sqrt_area_key,
    aggfunc='median',
    observed=True)
print('Median [Area]^(1/2):')
print(pivot)
print()

exm_median = pivot[expanded_cat_name].item()
unexm_median = pivot[unexpanded_cat_name].item()
ex_factor = exm_median / unexm_median

count_pivot = nuc_size_df.pivot_table(
    columns=category_key,
    values=sqrt_area_key,
    aggfunc='count',
    observed=True)

```

```
count_pivot = count_pivot.rename(index={sqrt_area_key: 'n'})
print('Num Data Points:')
print(count_pivot)
print()
```

```
exm_n = count_pivot[expanded_cat_name].item()
unexm_n = count_pivot[unexpanded_cat_name].item()
print(f'Expansion Factor: {ex_factor:.1f}x')
```

```
Median [Area]^(1/2):
category      expanded      normal
[Area]^(1/2)  34.59245   8.826096
```

```
Num Data Points:
category  expanded  normal
n          54      148
```

Expansion Factor: 3.9x

### 3.7 Create and Save Figure

```
fig = plt.figure(
    figsize=(2.5, 4.5),
    dpi=300)

palette_p = sns.color_palette(
    [sns.color_palette('pastel')[0],
     sns.color_palette('pastel')[3]])
palette_n = sns.color_palette(
    [sns.color_palette()[0],
     sns.color_palette()[3]])

ax = sns.boxplot(
    data=nuc_size_df,
    x=category_key,
    y=sqrt_area_key,
    hue=category_key,
    palette=palette_p,
    order=[unexpanded_cat_name, expanded_cat_name],
    hue_order=[unexpanded_cat_name, expanded_cat_name],
    linecolor='k',
```

```

    medianprops={'lw': 1.5})
sns.swarmplot(
    data=nuc_size_df,
    x=category_key,
    y=sqrt_area_key,
    hue=category_key,
    size=2.5,
    alpha=0.3,
    ax=ax,
    palette=palette_n,
    color='k',
    order=[unexpanded_cat_name, expanded_cat_name],
    hue_order=[unexpanded_cat_name, expanded_cat_name])
ax.axhline(
    y=unexm_median,
    color=palette_n[0],
    ls='--',
    lw=1,
    alpha=0.75
)
ax.axhline(
    y=exm_median,
    color=palette_n[1],
    ls='--',
    lw=1,
    alpha=0.75
)
sns.despine(right=True, bottom=True, top=True)
ax.set_ylim(0, 50)
ax.spines['bottom'].set_bounds(0, 1)
ax.set_yticks(
    [0,
     unexm_median,
     25,
     exm_median,
     50])
ax.set_yticklabels(
    ['0',
     f'{unexm_median:.1f}',
     '25',
     f'{exm_median:.1f}',

```

```

    '50'])
ax.set_ylabel('\u221A nucleus area (\u00b5m)')
ax.set_xlabel(None)
ax.annotate(
    f'{ex_factor:.1f}\u00D7',
    xy=(0.5, exm_median),
    xytext=(0.5, 20.0),
    arrowprops=dict(
        arrowstyle='->',
        color='black',
    ),
    annotation_clip=False,
    ha='center',
    va='center',
    size=14,
    weight='bold',
)
ax.annotate(
    f'{ex_factor:.1f}\u00D7',
    xy=(0.5, unexm_median),
    xytext=(0.5, 20.0),
    arrowprops=dict(
        arrowstyle='-',
        color='black',
    ),
    annotation_clip=False,
    ha='center',
    va='center',
    alpha=0,
    size=14,
    weight='bold',
)
ax.text(
    0, 2,
    f'n = {unexm_n}',
    ha='center',
    va='baseline',
    size=7.5,
    style='oblique',
    alpha=0.75)
ax.text(

```

```

1, 2,
f'n = {exm_n}',
ha='center',
va='baseline',
size=7.5,
style='oblique',
alpha=0.75)
ax.tick_params(axis='x', length=0, pad=-4)

x = np.array([31, 250, 500, 4000])

def save_figures(filename=None, figs=None, dpi=200, fmt='pdf', directory=None,
                 add_filename_timestamp=True, stamp_kwargs=None):
    from matplotlib.backends.backend_pdf import PdfPages

    if filename is None:
        filename = 'all_figures'

    if add_filename_timestamp:
        filename = filename + '_' + time.strftime('%y%m%d-%H%M')

    if figs is None:
        figs = [plt.figure(n) for n in plt.get_fignums()]
    else:
        try:
            _ = iter(figs)
        except TypeError:
            figs = [figs, ]

    if stamp_kwargs:
        [stamp_fig(fig, **stamp_kwargs) for fig in figs]

    directory = 'figures' if directory is None else directory
    try:
        os.mkdir(directory)
    except FileExistsError:
        pass

    file_path = os.path.join(directory, filename)
    if fmt == 'pdf':
        file_path += '.pdf'

```




```

with PdfPages(file_path) as pp:
    for fig in figs:
        fig.savefig(pp, format='pdf', dpi=dpi)
else:
    for i, fig in enumerate(figs):
        fig.savefig('%s_%d.%s' % (file_path, i, fmt), format=fmt, dpi=dpi)

plt.tight_layout()
save_figures(
    filename=figure_name,
    add_filename_timestamp=False,
    directory=figure_dir)

```



## 4 Notebook Conversion

This notebook can be converted to a .pdf via the following **pandoc** (<https://pandoc.org/>) command:

*for Unix-style command line:*

```
pandoc -o "sf3.pdf" --pdf-engine=xelatex --highlight-style tango \  
-V geometry:landscape \  
-V geometry:top=0.75in,bottom=0.75in,left=1in,right=2in \  
supplementary_file_01-nuclei-area-quant.ipynb
```

*for Windows Command Prompt:*

```
pandoc -o "sf3.pdf" --pdf-engine=xelatex --highlight-style tango ^  
-V geometry:landscape ^  
-V geometry:top=0.75in,bottom=0.75in,left=1in,right=2in ^  
supplementary_file_01-nuclei-area-quant.ipynb
```